




CT26/Human SSTR2 Stable Cell Line

Product Code	CSB-SC022725HU
Target Name	Adherent
Storage	Frozen in liquid nitrogen or stored at -80°C
Uniprot No.	P30874
Tested Applications	Binding assay by FACS and cell based ELISA.
Storage Buffer	35% FBS + 55% RPMI 1640 + 10% DMSO
Resistance	Puromycin
Intended Use	RPMI 1640 + 10% FBS + 10μg/mL Puromycin
Tag Info	Tag-Free

Image



Untransfected CT26 cells (green line) and transfected Human SSTR2 CT26 stable cells (red line) were stained with anti-SSTR2 antibody (CSB-RA022725MA01HU) (2μg/1*10⁶cells), washed and then followed by FITC-conjugated anti-Human IgG Fc antibody and analyzed with flow cytometry.