




CHOK1/Human CD93 Stable Cell Line

Product Code	CSB-SC865099HU
Target Name	Adherent
Storage	Frozen in liquid nitrogen or stored at -80°C
Uniprot No.	Q9NPY3
Tested Applications	Binding assay by FACS and cell based ELISA.
Storage Buffer	35% FBS + 55% F-12K + 10% DMSO
Resistance	Puromycin
Intended Use	F-12K + 10% FBS + 10µg/mL Puromycin
Tag Info	Tag-Free

Image



Untransfected CHO-K1 cells (green line) and transfected Human CD93 CHOK1 Stable cells (red line) were stained with anti-CD93 antibody(CSB-RA865099MA1HU) (2µg/1*10⁶cells), washed and then followed by FITC-conjugated anti-Human IgG Fc antibody and analyzed with flow cytometry.