





HEK293T/Human SARM1 Stable Cell Line

Product Code	CSB-SC750971HU
Target Name	Adherent
Storage	Frozen in liquid nitrogen or stored at -80°C
Uniprot No.	Q6SZW1
Tested Applications	Binding assay by FACS and Western blot
Storage Buffer	35% FBS + 55% DMEM + 10% DMSO
Resistance	Puromycin
Intended Use	DMEM + 10% FBS + 2μg/mL Puromycin
Tag Info	C-terminal Flag-tagged

Image



Untransfected HEK293T cells (green line) and transfected Human SARM1 HEK293T Stable cells (red line) were stained with anti-Flag antibody (2μg/1*10⁶cells), washed and then followed by FITC-conjugated anti-Mouse IgG Fc antibody and analyzed with flow cytometry.



Western Blot

Positive WB detected in: Lane1?20 μg HEK293T cells, Lane2?40 μg Human SARM1 HEK293T Stable cells, Lane3?20 μg Human SARM1 HEK293T Stable cells

All lanes: Anti-Flag antibody at 1:1000

Secondary

Goat polyclonal to mouse IgG at 1/50000 dilution

Predicted band size: 77.5 KDa

Observed band size: 77 KDa