




# HEK293T/Human CDH17 Stable Cell Line

<b>Product Code</b>	CSB-SC613267HU2
<b>Target Name</b>	Adherent
<b>Storage</b>	Frozen in liquid nitrogen or stored at -80°C
<b>Uniprot No.</b>	Q12864
<b>Tested Applications</b>	Binding assay by FACS and cell based ELISA.
<b>Storage Buffer</b>	35% FBS + 55% DMEM + 10% DMSO
<b>Resistance</b>	Puromycin
<b>Intended Use</b>	DMEM + 10% FBS + 2µg/mL Puromycin
<b>Tag Info</b>	Tag-Free

## Image



Untransfected HEK293T cells (green line) and transfected Human CDH17 HEK293T stable cells (red line) were stained with anti-CDH17 antibody (CSB-RA613267MA1HU) (2µg/1\*10<sup>6</sup>cells), washed and then followed by FITC-conjugated anti-Human IgG Fc antibody and analyzed with flow cytometry.