





# Recombinant *Macaca fascicularis* Membrane spanning 4-domains A1 (MS4A1)-VLPs (Active)

<b>Product Code</b>	CSB-MP4516MOV
<b>Abbreviation</b>	Recombinant <i>Cynomolgus</i> monkey MS4A1 protein-VLPs (Active)
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	M4ZH26
<b>Form</b>	Lyophilized powder
<b>Storage Buffer</b>	Lyophilized from a 0.2 μm filtered PBS, 6% Trehalose, pH 7.4
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	<i>Macaca fascicularis</i> (Crab-eating macaque) ( <i>Cynomolgus</i> monkey)
<b>Biological Activity</b>	Measured by its binding ability in a functional ELISA. Immobilized <i>Macaca fascicularis</i> CD20 at 10 μg/mL can bind Anti-CD20 recombinant antibody (CSB-RA015007MA3HU), the EC50 is 10.65-15.26 ng/mL. The VLPs (CSB-MP3838) is negative control.
<b>Sequence</b>	MTTPRNSVNGTFPAEPMKGP IAMQPGPKLLRRMSSLVGPTQSFFMRESKAL GAVQIMNGLFHIALGGLLMIPAGIYAPICVTVWYPLWGGIMYIISGSLLAATEKNS RKCLVKGKMIMNSLSLFAAISGMILSIMDILNIKISHFLKMESLNFIRVHTPYINIYN CEPANPSEKNPSTQYCYSIQSLFLGILSVMLIFAFFQELVIAGIVENEWRTCS RPKSSVLLSAEEKKEQVIEIKEEVVGLTETSSQPKNEEDIEIPIQEEEEETET NFPEPPQDQESSPIENDSSP
<b>Source</b>	Mammalian cell
<b>Target Names</b>	MS4A1
<b>Expression Region</b>	1-297aa
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	C-terminal 10xHis-tagged (This tag can be tested only under denaturing conditions)
<b>Mol. Weight</b>	34.4 kDa
<b>Protein Length</b>	Full Length

**Image**



CSB-MP4516MOV is detected by Mouse anti-6\*His monoclonal antibody.



#### Activity

Measured by its binding ability in a functional ELISA. Immobilized *Macaca fascicularis* CD20 at 10 µg/mL can bind Anti-CD20 recombinant antibody (CSB-RA015007MA3HU), the EC<sub>50</sub> is 10.65-15.26 ng/mL. The VLPs (CSB-MP3838) is negative control.

#### Endotoxin

Less than 1.0 EU/ug as determined by LAL method.

#### Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

#### Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.