




Recombinant Epstein-Barr virus Protein LF2 (LF2)

Product Code	CSB-EP314485EFA
Abbreviation	Recombinant Epstein-Barr virus LF2 protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P0C725
Form	Liquid or Lyophilized powder
Storage Buffer	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose.
Product Type	Recombinant Epstein-Barr virus Protein LF2(LF2)
Immunogen Species	Epstein-Barr virus (strain B95-8) (HHV-4) (Human herpesvirus 4)
Purity	≥ 90% as determined by SDS-PAGE.
Sequence	MAEAYPGGAHAALASRRSSFRNSLRRLRPTEKPDTSFMRGVWKYEIFPSYVR VTNKQVLQLDAQCQELPPCPSVGQILSFKLPSFSFNTTTYGSRYFTVAFLFFGA EDNEVFLKPFFVMHSDQDIVLSVLNPRSLFIEKGKFTWYIVPIRLVKNPYLYLQIL PGQSDIQLTRSCTQSGDKLNTSEPQIFLSGSPVTSQDECLPYLLAQHTPPFLKS YARIHTFPGKVCPVNAIRRGKGYVRVSVDPDLKREGPLNVKVGMTLLDDVIAA FRYNPYPKSHWRWDGESTDIRYFGSPVIIPPNFITELEYNNTYEAPLSSKITAVV VSHSSNPVFYVYPQEWKPGQTLKLTVRNISNNPITVTGQSMAQAFFIYAGDPS ISTIMRRYIQRQGCALTLPGNIVVESSSLPTFERINKTFNGNIVASEGTL
Research Area	Others
Source	E.coli
Target Names	LF2
Expression Region	1-429aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 10xHis-tagged and C-terminal Myc-tagged
Mol. Weight	53.5 kDa
Protein Length	Full Length
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.