



# Recombinant Bordetella pertussis BrkA autotransporter (brkA)

<b>Product Code</b>	CSB-EP673803BUA(A5)
<b>Abbreviation</b>	Recombinant Bordetella pertussis brkA protein
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q45340
<b>Form</b>	Liquid or Lyophilized powder
<b>Storage Buffer</b>	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose.
<b>Product Type</b>	Recombinant Bordetella pertussis BrkA autotransporter(brkA)
<b>Immunogen Species</b>	Bordetella pertussis (strain Tohama I / ATCC BAA-589 / NCTC 13251)
<b>Purity</b>	Greater than 85% as determined by SDS-PAGE.
<b>Sequence</b>	MYLDRFRQCPSSLQIPRSAWRLHALAAALALAGMARLAPAAAQAPQPPVAGA PHAQDAGQEGEFDHRDNTLIAVFDDGVGINLDDDPDELGETAPPTLKDIHISVE HKNPMSKPAIGVRVSGAGRALTLGASTIDATEGGIPAVVRRGGTLELDGVTVA GGEGMEPMTVSDAGSRLSVRGGVLGGEAPGVGLVRAAQGGQASIIDATLQSI LGPALIADGGSISVAGGSIDMDMGPFPPLPGAPLAAHPPLDRVA AVHA GQDGKVTLREVALRAHGPPQATGVYAYMPGSEITLQGGTVSVQGGDDGAGVVA GAGLLDALPPGGTVRLDGTTVSTDGANTDAVLVRGDAARAEVVNTVLR TAKSL AAGVSAQHGGRTLRQTRJETAGAGAEGISVLGFEPQSGSGPASVDMQGGSI TTTGNRAAGIALTHGSARLEGVAVRAEGSGSSAAQLANGTLVVSAGSLASAQS GAISVTDTPKLMPGALASSTVSVRLTDGATAQGGNGVFLQQHSTIPVAVALE SGALARGDIVADGNKPLDAGISLSVASGAAWHGATQVLQSATLGKGGTWVVN ADSRVQDMSMRGGRVEFQAPAPEASYKTLTLQTLTDGNGVFLNTNVAAGQN DQLRVTGRADGQHRVLRNAGGEADSRGARLGLVHTQGQGNATFRLANVGK AVDLGTWRYSLAEDPKTHVWSLQRAGQALSGAANA AVNAADLSSIALAESNA LDKRLGELRLRADAGGPWARTF SERQQISNRHARAYDQTVSGLEIGLDRGWS ASGGRWYAGLLGYTYADRTYPGDGGGKVKGLHVGGYAAAYVGDGGYYLDT VLRLLGRYDQQYNIAGTDGGRVTADYRTSGAAWSLEGGRRFELPNDWFAEPQ AEVMLWRTSGKRYRASNGLRVKVDANTATLGRGLRFRRIALAGGNIVQPY ARLGWTQEFKSTGDVRTNGIGHAGAGRHRVRELGAGVDAALGKGHNL YASY EYAAGDRINIPWSFHAGYRYSF
<b>Research Area</b>	Others
<b>Source</b>	E.coli
<b>Target Names</b>	brkA
<b>Expression Region</b>	1-1010aa




**Notes** Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

**Tag Info** N-terminal 10xHis-tagged and C-terminal Myc-tagged

**Mol. Weight** 110.8 kDa

**Protein Length** Full Length

**Image**



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

**Reconstitution** We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

**Shelf Life** The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.