




Recombinant Lordsdale virus Capsid protein (ORF2), Biotinylated

Product Code	CSB-MP346571LIV(A4)-B
Abbreviation	Recombinant Lordsdale virus ORF2 protein, Biotinylated
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P54635
Form	Liquid or Lyophilized powder
Storage Buffer	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose.
Product Type	Recombinant Lordsdale virus Capsid protein (ORF2), Biotinylated
Immunogen Species	Lordsdale virus (strain GII/Human/United Kingdom/Lordsdale/1993) (Human enteric calicivirus) (Hu/NV/LD/1993/UK)
Purity	≥ 90% as determined by SDS-PAGE.
Sequence	MKMASNDANPSDGSAAANLVPEVNNEVMALEPVVGAAIAAPVAGQQNVIDPWI RNNFVQAPGGEFTVSPRNAPGEILWSAPLGPDLNPYLSHLSRMYNGYAGGFE VQVILAGNAFTAGKVIFAAVPPNFPTEGLSPSQVTMFPHIIVDVRQLEPVLIP DVRNRFYHYNQANDSTLKIAMLTYPLRANNAGDDVFTVSCRVLTRPSPDFDF IFLVPPTVESRTKPFTVPVLTVEEMSNSRFPIPLEKLYTGPSSAFVVQPQNGRC TTDGVLLGTTQLSAVNICNFRGDVTHIAGSHDYTMNLASQNWSNYDPTEEIPA PLGTPDFVGKIQGLLTQTTRADGSTRAHKATVSTGSVHFTPKLGSVQFTTDTN NDFQAGQNTKFTPVGVIQDGDH HQNEPQQWSLPNYSGRTGHNVHLAPAVAP TFPGEQLLFFRSTMPGCSGYPNMNLDCLLPQEWVLHFYQEAPPAQSDVALLR FVNPD TGRVLFECKLHKSGYITVAHTGPYDLVLPNGYFRFDSWVNQFYTLAP MGNGTGRRRAL
Research Area	Others
Source	Mammalian cell
Target Names	ORF2
Expression Region	1-539aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 10xHis-tagged and C-terminal Myc-avi-tagged
Mol. Weight	64.3 kDa
Protein Length	Full Length



Image



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.