



Recombinant Human helicase MOV-10 (MOV10)

Product Code	CSB-EP862068HU(A4)
Abbreviation	Recombinant Human MOV10 protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9HCE1
Form	Liquid or Lyophilized powder
Storage Buffer	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose.
Product Type	Recombinant Human helicase MOV-10(MOV10)
Immunogen Species	Homo sapiens (Human)
Purity	≥ 90% as determined by SDS-PAGE.
Sequence	MPSKFSCRQLREAGQCFESFLVVRGLDMETDRERLRTIYNRDFKISFGTPAPG FSSMLYGMKIANLAYVTKTRVRRFFRLDRWADVRFPEKRRMKLGS DISKHHKSL LAKIFYDRAEYLHGKHGVDVEVQGPHEARDGQLLIRLDLNRKEVLT LRLRNGG TQSVTLTHLFP LCRTPQFAFY NEDQELPCPLGPGECYELHVHCKTSFVGYFPA TVLWELLGPGESGSEGAGTFYIARFLAAVAHSPLAAQLKPMTPFKRTRITGNP VVTNRIEEGERPDRAKGYDLELSMALGTY YPPPRLRQLL PMLLQGTSIFTAPKE IAEIKAQLETALKWRNYEVKLRLLLHLEELQMEHDIRHYDLESVPMTWDPVDQ NPRLLTLEVPGVTESRPSVLRGDHLFALLSSETHQEDPITYKGFVHKVELDRVK LSFSMSLLSRFVDGLTFKVNFTFN RQPLRVQHRAELTGRWLLWPMLFPVAP RDVPLLPDVK LKLYDRSLESNPEQLQAMRHIVTGTTRPAPYIIFGPPGTGKT V TLVEAIKQVVKHLPKAHILACAPSNSGADLLCQRLRVHLPSSIYRLLAPSRDIRM VPEDIKPCCNWDAKKGEYVFP AKKKLQEYRVLITTLITAGRLVSAQFPIDHFTHI FIDEAGHCMEPESLVAIAGLMEVKETGDPGGQLVLAGDPRQLGPVLRSP L TQK HGLGYSLLERLLTYNSLYKKGPDGYDPQFITKLLRNYRSHPTILDIPNQLYYEGE LQACADVDRERFCRWAGLPRQGFPIIFHGVMGKDEREGNSPSFFNP EEAAT VTSYLKLLLAPSSKKGKARLS PRSVGVISPYRKQVEKIRYCITKLDREL RGLDDI KDLKVGSV EEFQGGQERSVILISTVRSSQS FVQLDLDFNLGFLKNPKRFNVA VTR AKALLIIVGNPLLLGHDPDWKVFLEFCKENGGYTGC PFP AKLDLQGGQNLLQG LSKLS PSTSGPHSHDYLPQEREGEGGLSLQVEPEWRNEL
Research Area	Epigenetics and Nuclear Signaling
Source	E.coli
Target Names	MOV10
Expression Region	1-1003aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.




Tag Info N-terminal 6xHis-tagged

Mol. Weight 119.6 kDa

Protein Length Full Length

Image



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.