




Recombinant Enterobacteria phage T4 Single-stranded DNA-binding protein (32)

Product Code	CSB-EP022706EDZe1
Abbreviation	Recombinant Enterobacteria phage T4 Single-stranded DNA-binding protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P03695
Form	Liquid or Lyophilized powder
Storage Buffer	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose.
Product Type	Recombinant Enterobacteria phage T4 Single-stranded DNA-binding protein(32)
Immunogen Species	Enterobacteria phage T4 (Bacteriophage T4)
Purity	Greater than 90% as determined by SDS-PAGE.
Sequence	MFKRKSTAELAAQMAKLNKNGKGFSSSEDKGEWKLKLDNAGNGQAVIRFLPSKN DEQAPFAILVNHGFKKNGKWIETCSSTHGDYDSCPVCQYISKNDLYNTDNKE YSLVKRKTSYWANILVVKDPAAPENEGKVKYRFGKKIWDKINAMIAVDVEMG ETPVDVTCPWEGANFVLKVKQVSGFSNYDESKFLNQS AIPNIDDES FQKELFE QMVDLSEMTSKDKFKSFEELNTKFGQVMGTAVMGGAAATAAKKADKVADDL DAFNVDDFNTKTEDDFMSSSSGSSSSADDTDLDDLLNDL
Research Area	others
Source	E.coli
Target Names	32
Expression Region	1-301aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag-Free
Mol. Weight	33.5 kDa
Protein Length	Full Length
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.