




Recombinant Mumps virus Fusion glycoprotein F0 (F), partial

Product Code	CSB-EP322556MJM
Abbreviation	Recombinant Mumps virus Fusion glycoprotein F0, partial
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P19716
Form	Liquid or Lyophilized powder
Storage Buffer	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose.
Product Type	Recombinant Mumps virus Fusion glycoprotein F0(F),partial
Immunogen Species	Mumps virus (strain SBL-1) (MuV)
Purity	≥ 90% as determined by SDS-PAGE.
Sequence	FAGIAIGIAALGVATAAQVTA AVSLVQAQTNARAIAAMKNSIQATNRAVFEVKEG TQQLAIAVQAIQDHINTIMNTQLNNMSCQILDNQLATSLGLYLTELTTVFQPQLIN PALSPISIQALRSLLGSMTPAVVQATLSTSISAAEILSAGLMEGQIVSVLLDEMQ MIVKINVPTIVTQSNALVIDFYSSIFINNQESIIQLPDRILEIGNEQWRYPAKNCK LTRHHIFCQYNEAERLSLETKLCLAGNISACVFSSIAGSYMRRFVALDGTIVANC RSLTCLCKSPSYPIYQPDHHA VTTIDLTSCQTLSLDGLDFSIVSLSNITYAENLTI SLSQTINTQPIDISTELSKVNASLQNAVKYIKESNHQLQSVSVSSKIGAI
Research Area	Others
Source	E.coli
Target Names	F
Expression Region	103-486aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 10xHis-SUMO-tagged
Mol. Weight	56.0 kDa
Protein Length	Partial
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.