




Recombinant Pseudomonas putida RNA polymerase sigma factor RpoS (rpoS)

Product Code	CSB-EP4767GTH
Abbreviation	Recombinant Pseudomonas putida rpoS protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q88ME8
Form	Liquid or Lyophilized powder
Storage Buffer	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose.
Product Type	Recombinant Pseudomonas putida RNA polymerase sigma factor RpoS(rpoS)
Immunogen Species	Pseudomonas putida (strain ATCC 47054 / DSM 6125 / CFBP 8728 / NCIMB 11950 / KT2440)
Purity	≥ 90% as determined by SDS-PAGE.
Sequence	MALSKEVPEFDIDDDVLLMETGIVLETDVVSDEPAVSSVRTRAKSGSSLKQHK YIDYSRALDATQLYLNEIGFSPLLSPEEEVHFARLSQKGDPAGRKRMIESNLR LVVKIARRYVNRGLSLLDLIEEGNLGLIRAVEKFDPERGFRFSTYATWWIRQTIER AIMNQTRTIRLPIHVVKELNVYLRAARELTQKLDHEPSPEEIATLLEKPVAEVKR MLGLNERVSSVDVSLGPDSDKTLTDDRPTDPCCELLQDDDLQSIDSQWLG ELTDKQREVVVRRFGLRGHESSTLEDVGLIQLTRERVRQIQVEGLKRLREILE KNGLSSESLEFQ
Research Area	Others
Source	E.coli
Target Names	rpoS
Expression Region	1-335aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 10xHis-tagged
Mol. Weight	44.2 kDa
Protein Length	Full Length
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.