




Recombinant Vaccinia virus Intermediate transcription factor 3 large subunit (VITF3L)

Product Code	CSB-EP324999VAA
Abbreviation	Recombinant Vaccinia virus VITF3L protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P20998
Form	Liquid or Lyophilized powder
Storage Buffer	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose.
Product Type	Recombinant Vaccinia virus Intermediate transcription factor 3 large subunit(VITF3L)
Immunogen Species	Vaccinia virus (strain Copenhagen) (VACV)
Purity	≥ 90% as determined by SDS-PAGE.
Sequence	MDNLFTFLHEIEDRYARTIFNFHLISCDEIGDIYGLMKERISSEDMFDNIVYNKDI HHAIKKLVCYCDIQLTKHIINQNTYPVFNDSSQVKCCHYFDINSDNSNISSRTVEIF EREKSSLVSYIKTTNKKRKVNYGEIKKTVHGGTNANYFSGKKSDEYLSTTVRS NINQPWIKTISKMRVDIINHISIVTRGKSSILQTIEIIFTNRTC VKIFKDSTMHIILSK DKDEKGCIH MIDKLFYVYYNLFLLFEDI IQNEYFKEVANVNVNHLTATALDEKLFL IKKMAEHDVYGVS NF KIGMFNLTFIKSLDHTVFP SLLDEDSKIKFFKGKKLNI VAL RSLEDCIN YVT KSENMIEMMKERSTILNSIDIETESVDRLKELLLK
Research Area	Others
Source	E.coli
Target Names	VITF3L
Expression Region	1-382aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	C-terminal 6xHis-tagged
Mol. Weight	51.5 kDa
Protein Length	Full Length
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.