




# Recombinant Human Proprotein convertase subtilisin/kexin type 9 (PCSK9) (D374Y), partial

<b>Product Code</b>	CSB-MP017647HU(M)
<b>Abbreviation</b>	Recombinant Human PCSK9 protein, partial
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q8NBP7
<b>Form</b>	Liquid or Lyophilized powder
<b>Storage Buffer</b>	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose.
<b>Product Type</b>	Recombinant Human Proprotein convertase subtilisin/kexin type 9(PCSK9)(D374Y),partial
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	Greater than 90% as determined by SDS-PAGE.
<b>Sequence</b>	QEDEDGDYEELVLA LRSEEDGLAEAPEHGTTATFHRC AKDPWRLPGTYVVVL KEETHLSQSERTARRLQAQAARRGYLTKILHVFHGLLPGLVKMSGDLLELALK LPHVDYIEEDSSVFAQSIPWNLERITPPRYRADEYQPPDGGSLVEVYLLDTSIQ SDHREIEGRVMVTD FENVPEEDGTRFHRQASKCD SHGTHLAGVVSGRDAGV AKGASMRSLRVLN CQGKGT VSGTLIGLEFIRKSQLVQPVG PLVLLPLAGGYS RVLNAACQRLARAGVVLVTAAGNFRDDACLYSPASAPEVITVGATNAQDQPVT LGT LGTNFGRCVDL FAPGEDIIGASSY CSTCFVSQSGTSQAAAHVAGIAAMML SAPELTLAELRQRLIHFS AKDVINEAWFPEDQ RVLTPNLVAALPPSTHGAGW QLFCRTVWSAHSGPTRMATAVARCAPDEELLSCSSF SRSGKRRGERMEAQG GKLV CRAHNAFGGEGVYAIARCCLLPQANCSVHTAPPAEASMGTRVHCHQQ GHVLTGCSSHWEVD LGTHKPPVLRPRGQPNQC VGHREASIHASCCHAPGL ECKVKEHGIPAPQE QTVACEEGWTLTGCSALPGTSHVLGAYAVDNTCVVRS RDVSTTGSTSEGAVTAVAICCRSRHLAQASQELQ
<b>Research Area</b>	Cardiovascular
<b>Source</b>	Mammalian cell
<b>Target Names</b>	PCSK9
<b>Expression Region</b>	31-692aa
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	C-terminal 6xHis-tagged
<b>Mol. Weight</b>	73.3 kDa

**Protein Length**

Partial

**Image**

(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

**Reconstitution**

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

**Shelf Life**

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.