




Recombinant Vaccinia virus Soluble interferon gamma receptor OPG193 (OPG193)

Product Code	CSB-MP322719VAAAd9
Abbreviation	Recombinant Vaccinia virus Soluble interferon gamma receptor OPG193 protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P21004
Form	Liquid or Lyophilized powder
Storage Buffer	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose.
Product Type	Recombinant Vaccinia virus Soluble interferon gamma receptor B8(B8R)
Immunogen Species	Vaccinia virus (strain Copenhagen) (VACV)
Purity	Greater than 90% as determined by SDS-PAGE.
Sequence	SIHAKITSYKFEVSNFDSKIEWTGDGLYNISLKNYGIKTWQTMVYTNVPEGTYDIS AFPKNDFVSWVKFEQGDYKVEEYCTGLCVEVKIGPPTVTLTEYDDHINLYIEH PYATRGSKKIPIYKRGDMCDIYLLYTANFTFGDSEEPVTYDIDDYDCTSTGCSID FATTEKVCVTAQGATEGFLEKITPWSSEVCLTPKKNVYTCAIRSKEDVPNFKDK MARVIKRKFNKQSQSYLTKFLGSTSNDVTTFLSMLNLTKYS
Research Area	Others
Source	Mammalian cell
Target Names	OPG193
Expression Region	14-272aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	C-terminal hFc1-tagged
Mol. Weight	58.5 kDa
Protein Length	Full Length of Mature Protein
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.