




ERBB2 Recombinant Monoclonal Antibody

Product Code	CSB-RA007763MA1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P04626
Immunogen	Recombinant Human ERBB2 protein
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	Affinity-chromatography
Isotype	hIgG1
Clonality	Monoclonal
Product Type	Recombinant Antibody
Immunogen Species	Homo sapiens (Human)
Research Area	Cancer;Tags & Cell Markers?Immunology;Signal transduction
Target Names	ERBB2
Clone No.	6B3

Image



The Binding Activity of Human ERBB2 with Anti-ERBB2 Recombinant Antibody
 Activity: Measured by its binding ability in a functional ELISA. Immobilized Human ERBB2 (CSB-MP007763HU) at 2 µg/mL can bind Anti-ERBB2 recombinant antibody, the EC₅₀ is 2.455-3.565 ng/mL.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.