



Recombinant *Macaca fascicularis* Cadherin 17 (CDH17), partial (Active)


Product Code	CSB-MP4664MOV
Abbreviation	Recombinant <i>Cynomolgus</i> monkey CDH17 protein, partial (Active)
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	A0A2K5X8I8
Form	Lyophilized powder
Storage Buffer	Lyophilized from a 0.2 µm sterile filtered PBS, 6% Trehalose, pH 7.4
Product Type	Recombinant Protein
Immunogen Species	<i>Macaca fascicularis</i> (Crab-eating macaque) (<i>Cynomolgus</i> monkey)
Biological Activity	Measured by its binding ability in a functional ELISA. Immobilized <i>Macaca fascicularis</i> CDH17 at 2 µg/mL can bind Anti-CDH17 recombinant antibody (CSB-RA613267MA3HU). The EC50 is 2.182-2.529 ng/mL.
Purity	≥ 90% as determined by SDS-PAGE.
Sequence	KILSEIGDLLQIPQLCFFPYSLICFFPQFKANPPAVTFELTGETDNIFKIEQEGLL YYTKALDRETRSTHNLQVAALDANGAIVEGVPITIEVKDVNDNRPTFLQSKYE GSVRQNSRPGKPFLLYVNATDLDDPATPNGQLSYQIVIQLPMINNVMYFQINNK TGGISLTREGSQELNPAKNPSYNLVISVKDMGGQSENSFSDTTSVDIIVTENIW KAPEPVEMVENSTDPHPKITQVRWNDPGAQYSLVDKEKLPFRPFPSIDQEGDIY VTQPLDREEKDAYVFYAVAKDEYKPLSYPLEIHVKVQDINDNPPTCPSPTVF EVQENERLGNISGALTAHDSDEENTANLLNYRIVEQTPKLPMDGLFLIQTYAG MLQLAKQSLKKQDTPQYNLTIEVSDKDFKTLCFVQINVIDINDQIPIFEKSDYGNL TLAEDTNVGSTILTIQATDADEPFTGSSKILYHVIKGDSEGRGLGVDTDPHTNTGY VIIKKPLDFETAISNIVFKAENPEPLVFGVTYNASSFAKFTLVTDVNEAPEFSQ YVFQAKVSEDVAIGTKVGNVTAKDPEGLDISYSLRGDTRGWLKIDHVTGEIFSV APLDREAGSPYRVQVVATEVGGSSLSVVSQFHLILTDVNDNPPRLAKDYMDLY FCHPLSAPGSLIFEATDDDQHLFRGPHFTFSIASESLQNDWQVSKINGTHARLS TRHTDFEEKEYVVSIRINDGGRPPLESTVSLTVTFCSCGEDGCFRPAHQPGI PT
Source	Mammalian cell
Target Names	CDH17
Expression Region	30-792aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	C-terminal 10xHis-tagged




Mol. Weight 86.4 kDa

Protein Length Partial

Image



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.



Activity
Measured by its binding ability in a functional ELISA. Immobilized *Macaca fascicularis* CDH17 at 2 µg/ml can bind Anti-CDH17 recombinant antibody (CSB-RA613267MA3HU). The EC₅₀ is 2.182-2.529 ng/mL.

Endotoxin Less than 1.0 EU/ug as determined by LAL method.

Reconstitution We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.