





# Recombinant *Macaca fascicularis* zymogen granule protein 16 homolog B (ZG16B) (Active)

<b>Product Code</b>	CSB-MP4727MOV
<b>Abbreviation</b>	Recombinant Cynomolgus monkey ZG16B protein (Active)
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	G7Q0A1
<b>Form</b>	Lyophilized powder
<b>Storage Buffer</b>	Lyophilized from a 0.2 µm filtered 20 mM Tris-HCl, 0.5 M NaCl, 6% Trehalose, pH 8.0
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	<i>Macaca fascicularis</i> (Crab-eating macaque) (Cynomolgus monkey)
<b>Biological Activity</b>	①Measured by its binding ability in a functional ELISA. Immobilized <i>Macaca fascicularis</i> ZG16B at 2 µg/mL can bind Anti-ZG16B recombinant antibody (CSB-RA836195MA1HU), the EC50 is 20.15-28.98 ng/mL.
<b>Purity</b>	Greater than 90% as determined by SDS-PAGE.
<b>Sequence</b>	GKMYGPGGGKYFTTTEDYDHEITGLRVSVGLLLVKSQVVKLGDTWDVKQGAS GGNTQEVTLPGEYITKVFVAFQTFLRGMVLYTSKDRFTFYFGKLDGQIFSVYP SQEGQVLVGIYGQYGLLGIKISIGFEWNYPLEEPTTEPPVTVT
<b>Research Area</b>	Cancer
<b>Source</b>	Mammalian cell
<b>Target Names</b>	ZG16B
<b>Expression Region</b>	23-169aa
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	C-terminal 10xHis-tagged
<b>Mol. Weight</b>	17.6 kDa
<b>Protein Length</b>	Full Length of Mature Protein
<b>Image</b>	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.



**Activity**  
 Measured by its binding ability in a functional ELISA. Immobilized *Macaca fascicularis* ZG16B at 2 $\mu$ g/mL can bind Anti-ZG16B recombinant antibody (CSB-RA836195MA1HU)? the EC<sub>50</sub> is 20.15-28.98 ng/mL.

**Endotoxin**

Less than 1.0 EU/ $\mu$ g as determined by LAL method.

**Reconstitution**

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

**Shelf Life**

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.