




Recombinant Mouse Putative ATP-dependent RNA helicase DHX33 (Dhx33)

Product Code	CSB-EP767332MO
Abbreviation	Recombinant Mouse Dhx33 protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q80VY9
Form	Liquid or Lyophilized powder
Storage Buffer	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose.
Product Type	Recombinant Mouse Putative ATP-dependent RNA helicase DHX33(Dhx33)
Immunogen Species	Mus musculus (Mouse)
Sensitivity	Not Test
Purity	Greater than 85% as determined by SDS-PAGE.
Sequence	MPEEASLPPAKRFRPGSCPPGRRVVMLLTAGGGGGAGGGRRQTPPLAQPSA SPYREALELQRRSLPIFRARGQLLAQLRNLDAVLIGETGSGKTTQIPQYLYEG GISRQGIIAVTQPRRVA AISLATRV SDEKRT ELGKLVGYTVRFEDVTS EDTRIKF LTDGMLLRE AISDSLLRKYSCVILDEAHERTIHTDVLFGVVKTAQKRRKELGKLP LKVIVMSATMDVDLFSQYFNRA PVLYLEGRQHPIQIFYTKQPQQDY LHAALVSV FQIHQEAPASQDILVFLTGQEEIEAMSKTCRDIARHLPDGCPSMLVLPYASLP YSQQLRVFQGAPKGYRKVIISTNIAETSITITGIKYVVDTGMMVKAKKYNPDSGLE VLAVQRVSKTQAWQRTGRAGREDSGICYRLYTEDEFEFKFKMTVPEIQRCNL ASVILQLLAMKVPNVLTFDFMSKPSPDHIEAAIAQLDLLGALEHKDDQLTLTPIG RKMAAFPLEPRFAKTILLSSKFHCTEEILTIVSLLSVDSVLYNPPARRDEVQSVR KKFISSEGDHITLLNIYRTFKNIGGNKDWCKENFVNSKNMMLLVAEVRAQLREICL KMSMPIMSSRGDMESVRRRCMAHSLFMNTAELQTDGTYATTDTHTQPVAIHPS VLFHCKPACVVYTSLLYTNKCYMRDLCVVD AEWLYEAAPDYFRRKLRTARN
Research Area	Cell Biology
Source	E.coli
Target Names	Dhx33
Expression Region	1-698aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 6xHis-tagged
Mol. Weight	82.4 kDa

**Protein Length**

Full Length

Image

(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.