




Recombinant Human Hypoxia up-regulated protein 1 (HYOU1), partial

Product Code	CSB-EP896895HU
Abbreviation	Recombinant Human HYOU1 protein, partial
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9Y4L1
Form	Liquid or Lyophilized powder
Storage Buffer	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose.
Product Type	Recombinant Human Hypoxia up-regulated protein 1(HYOU1),partial
Immunogen Species	Homo sapiens (Human)
Sensitivity	Not Test
Purity	Greater than 90% as determined by SDS-PAGE.
Sequence	MVEEIGVELVVLDPDLPEDKLAQSVQKLQDLTLRDLEKQEREKAANSLEAFIF ETQDKLYQPEYQEVSTEEQREEISGKLSAASTWLEDEGVGATTVMLKEKLAEL RKLCQGLFFRVEERKKWPERLSALDNLLNHSSMFLKGARLIPEMDQIFTEVEM TTLEKVINETWAWKNATLAEQAKLPATEKPVLLSKDIEAKMMALDREVQYLLNK AKFTKPRPRPKDKNGTRAEPPLNASASDQGEKVIIPAGQTEDAEPISEPEKVE TGSEPGDTEPLELGGPGAPEQKEQSTGQKRPL
Research Area	Cancer
Source	E.coli
Target Names	HYOU1
Expression Region	695-994aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 6xHis-tagged
Mol. Weight	37.8 kDa
Protein Length	Partial
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.