



Recombinant Lactobacillus delbrueckii subsp. bulgaricus Chaperonin GroEL (groEL), partial

Product Code	CSB-EP625206LAQ
Abbreviation	Recombinant Lactobacillus delbrueckii subsp. bulgaricus groEL protein, partial
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q1G937
Form	Liquid or Lyophilized powder
Storage Buffer	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol.If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose.
Product Type	Recombinant Lactobacillus delbrueckii subsp. bulgaricus 60 kDa chaperonin(groL),partial
Immunogen Species	Lactobacillus delbrueckii subsp. bulgaricus (strain ATCC 11842 / DSM 20081 / BCRC 10696 / JCM 1002 / NBRC 13953 / NCIMB 11778 / NCTC 12712 / WDCM 00102 / Lb 14)
Sensitivity	Not Test
Purity	≥ 90% as determined by SDS-PAGE.
Sequence	ETELSVVEGMQFDRGYLSQYMVTDNDKMEADLENPYILITDKKISNIQDILPML QEIVQQGRSLLIADDVTGEALPTLVLNKIRGTFNVVAVKAPGFGDRRKEQLADI AALTGGTVISEDLGLELKDTQLSQLGQARRVTITKDSTTIVDGSGAKEAIQERV DTIRKQIEDTSSDFDKKKLQERLAKLTGGV
Research Area	Others
Source	E.coli
Target Names	groEL
Expression Region	182-374aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 10xHis-tagged and C-terminal Myc-tagged
Mol. Weight	26.3 kDa
Protein Length	Partial
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.