




TFRC Recombinant Monoclonal Antibody

Product Code	CSB-RA023441MA1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P02786
Immunogen	Recombinant Human TFRC protein
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	Affinity-chromatography
Isotype	hIgG1
Clonality	Monoclonal
Product Type	Recombinant Antibody
Immunogen Species	Homo sapiens (Human)
Research Area	Cancer; Cardiovascular; Immunology; Metabolism; Stem cells
Target Names	TFRC
Clone No.	7H12

Image



The Binding Activity of Human TFRC with Anti-TFRC Recombinant Antibody
 Activity: Measured by its binding ability in a functional ELISA. Immobilized Human TFRC (CSB-MP3648HU) at 2 µg/mL can bind Anti-TFRC recombinant antibody, the EC₅₀ is 3.305-8.220 ng/mL.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.