



PNPLA3 Antibody

Product Code	CSB-PA154215
Storage	upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NST1
Immunogen	Synthesized peptide derived from human PNPLA3
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	WB, ELISA; WB:1:1000-1:2000, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	PNPLA3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.