



# Recombinant Mouse Cap-specific mRNA (nucleoside-2-O-)-methyltransferase 2 (Cmtr2)

<b>Product Code</b>	CSB-EP804417MO
<b>Abbreviation</b>	Recombinant Mouse Cmtr2 protein
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q8BWQ4
<b>Form</b>	Liquid or Lyophilized powder
<b>Storage Buffer</b>	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose.
<b>Product Type</b>	Recombinant Mouse Cap-specific mRNA(nucleoside-2-O-)-methyltransferase 2(Cmtr2)
<b>Immunogen Species</b>	Mus musculus (Mouse)
<b>Sensitivity</b>	Not Test
<b>Purity</b>	≥ 85% as determined by SDS-PAGE.
<b>Sequence</b>	MSKRRKLPARQPACLETFS PDV LNDVSE LFAKSF SYRKPLDNEWQLPAPTESF SCGHLEFRALLDLKNSLNEVKNLLSDKKLDEWHRHTAFTNKAGKIISHVKKAVN AELCTQAWCKFQEILCSFPLIPQEA FQSGRLNSLHLCEAPGAFIASLNHYLKSH RFPCEWSWVANS LN PYHEANDNLRMITDDRLMANTLHCWYFGPDNTGDIMTL KYLTGLQDFLSGMSP IHLVTADGSFDCQGNPGEQEALVSSLHYCEAVTALITL GDGGSFVLKMFTLFEHCSVNLMYLLNCSFDQVHVFKPATSKAGNSEVYVVCL RYKGREAVQPLL SRMVLNFGTEMTRKALFP HHVIPKSF LERHEECCTFFHRYQ LETISENIRLFESMGTGEQERLNNLRDCAVQYFMQKFQLKPLSRNHVLVKKSN IGCSMNTKWFGQRNKYFKTYNERKMMETLSWKDKVAKGYFN SWAE EHTVYH PGQNSLLEGTASSLEYQSWQVLEGGKLPKVKCSPFCDGEILKTLNEAIEKSLG EALSVDKAVSSKQQYRCCPVFSEESVLS ELLRLTKCLPDEQGAEPSPVKCLL VGSPA VCDLQMPAPLEIQLVESVELTAFSCSLLHDGDPAYQHLFLDCLLHSLR RLHRGDVMVLPILSCFTRFMAGLTFVLHGCFR FITSCPTSLEPLRTC AVL LCIG YQNL PDAVFQFLQNVHDLLSKLLHPSAPRQILQFLPMEALLQG TLLDFLWDLNA AIAKRHLHLIIQGERDQAIGSLEL
<b>Research Area</b>	Others
<b>Source</b>	E.coli
<b>Target Names</b>	Cmtr2
<b>Expression Region</b>	1-767aa
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.




**Tag Info** N-terminal 10xHis-tagged

**Mol. Weight** 90.0 kDa

**Protein Length** Full length

**Image**



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

**Reconstitution**

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

**Shelf Life**

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.