




cry1Ia Antibody

| | |
|----------------------------|--|
| Product Code | CSB-PA673817XA01BDC |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | Q45752 |
| Immunogen | Recombinant Bacillus thuringiensis subsp. kurstaki cry1Ia protein (97-180aa) |
| Species Reactivity | Bacillus thuringiensis subsp. kurstaki |
| Tested Applications | ELISA, WB |
| Form | Liquid |
| Conjugate | Non-conjugated |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4 |
| Purification Method | Affinity-chromatography |
| Isotype | IgG |
| Clonality | Polyclonal |
| Product Type | Polyclonal Antibody |
| Immunogen Species | Bacillus thuringiensis subsp. kurstaki |
| Target Names | cry1Ia |

Image



Western Blot

Positive WB detected in: recombinant protein

All lanes: cry1Ia Antibody at 1:1000

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 35 kDa

Observed band size: 35 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.