




# G Antibody

<b>Product Code</b>	CSB-PA4065XA01
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	B7UCZ5
<b>Immunogen</b>	Recombinant vesicular stomatitis Indiana virus rVSV-G/GFP G protein (31-232aa)
<b>Species Reactivity</b>	Recombinant vesicular stomatitis Indiana virus rVSV-G/GFP
<b>Tested Applications</b>	ELISA, WB
<b>Form</b>	Liquid
<b>Conjugate</b>	Non-conjugated
<b>Storage Buffer</b>	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
<b>Purification Method</b>	Affinity-chromatography
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Product Type</b>	Polyclonal Antibody
<b>Immunogen Species</b>	Recombinant vesicular stomatitis Indiana virus rVSV-G/GFP
<b>Target Names</b>	G

## Image



### Western Blot

Positive WB detected in: recombinant protein

All lanes: G Antibody at 1:1000

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 35 kDa

Observed band size: 35 kDa

## Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.