




Early E3 18.5 kDa glycoprotein Antibody

Product Code	CSB-PA361279XA01HIL
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P04494
Immunogen	Recombinant Human adenovirus C serotype 5 Early E3 18.5 kDa glycoprotein protein (18-124aa)
Species Reactivity	Human adenovirus C serotype 5
Tested Applications	ELISA, WB
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Affinity-chromatography
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Human adenovirus C serotype 5 (HAdV-5) (Human adenovirus 5)
Target Names	N/A

Image



Western Blot

Positive WB detected in: recombinant protein

All lanes: Early E3 18.5 kDa glycoprotein

Antibody at 1:1000

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 30 kDa

Observed band size: 30 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.