




ASPZODRAFT_1728931 Antibody

Product Code	CSB-PA4716XA01DLU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	A0A1L9SK73
Immunogen	Recombinant Penicilliosis zonata CBS 506.65 ASPZODRAFT_1728931 protein (302-356aa)
Species Reactivity	Penicilliosis zonata CBS 506.65
Tested Applications	ELISA, WB
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Affinity-chromatography
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Penicilliosis zonata CBS 506.65
Target Names	ASPZODRAFT_1728931

Image



Western Blot

Positive WB detected in: recombinant protein
All lanes: ASPZODRAFT_1728931 Antibody at 1:1000

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 24 kDa

Observed band size: 24 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.