




ASPZODRAFT_1728931 Antibody

| | |
|----------------------------|--|
| Product Code | CSB-PA4716XA11DLU |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | A0A1L9SK73 |
| Immunogen | Recombinant Penicilliosis zonata CBS 506.65 ASPZODRAFT_1728931 protein (173-274aa) |
| Species Reactivity | Penicilliosis zonata CBS 506.65 |
| Tested Applications | ELISA, WB |
| Form | Liquid |
| Conjugate | Non-conjugated |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4 |
| Purification Method | Affinity-chromatography |
| Isotype | IgG |
| Clonality | Polyclonal |
| Product Type | Polyclonal Antibody |
| Immunogen Species | Penicilliosis zonata CBS 506.65 |
| Target Names | ASPZODRAFT_1728931 |

Image



Western Blot

Positive WB detected in: recombinant protein
All lanes: ASPZODRAFT_1728931 Antibody at 1:1000

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 15 kDa

Observed band size: 15 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.