



TaV-CP Antibody

Product Code	CSB-PA4642XA11GSU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9YK87
Immunogen	Synthetic Peptide (139-152aa)
Species Reactivity	Thosea asigna virus
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Affinity-chromatography
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Thosea asigna virus
Target Names	TaV-CP
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.