




CTAG2 Antibody

Product Code	CSB-PA029753
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O75638
Immunogen	Synthesized peptide derived from part region of human protein
Raised In	Rabbit
Species Reactivity	Human, Mouse, Rat
Tested Applications	ELISA, WB; Recommended dilution: WB:1:500-1:2000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	CTAG2

Image



Western blot analysis of lysates from DU145 cells, primary antibody was diluted at 1:1000, 4° over night

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.