



Recombinant Human Delta-1-pyrroline-5-carboxylate synthase (ALDH18A1)

Product Code	CSB-EP001564HUe1
Abbreviation	Recombinant Human ALDH18A1 protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P54886
Form	Liquid or Lyophilized powder
Storage Buffer	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose.
Product Type	Recombinant Human Delta-1-pyrroline-5-carboxylate synthase(ALDH18A1)
Immunogen Species	Homo sapiens (Human)
Sensitivity	Not Test
Purity	≥ 85% as determined by SDS-PAGE.
Sequence	<p> MLSQVYRCGFQPFNQHLLPWVKCTTVFRSHCIQPSVIRHVRSWSNIPFITVPLS RTHGKSFAHRSELKHAKRIVVKLGSVVTRGDECGLALGRLASIVEQVSVLQN QGREMMLVTSGAVAFGKQRLRHEILLSQSVRQALHSGQNQLKEMAIPVLEAR ACAAAGQSGLMALYEAMFTQYSICAAQILVTNLDFHDEQKRRNLNGTLHELLR MNIVPIVNTNDVVPPAEPNSDLQGVNVISVKDNDSLAARLAVEMKTDLLIVLS DVEGLFDSPPGSDDAKLIDIFYPGDQQSVTFGTKSRVGMGGMEAKVKAALWA LQGGTSVVIANGTHPKVSGHVITDIVEGKKVGTFFSEVKPAGPTVEQQGEMAR SGGFMLATLEPEQRAEIIHHLADLLTDQRDEILLANKKDLEEAEGRLAAPLLKRL SLSTSKLNSLAIGLRQIAASSQDSVGRVLRRTRIAKNLELEQVTVPIGVLLVIFES RPDCLPQVAALAIASGNLLLLKGGKEAAHSNRILHLLTQEALSIHGKVEAVQLV NTREEVEDLCRLDKMIDLIIPRGSSQLVRDIQKAAKGIPVMGHSEGICHMYVDS EASVDKVTRLVVRDSKCEYPACNALETLLIHRDLLRTPLFDQIIDMLRVEQVKIH AGPKFASYLTFSPSEVKSLRTEYGDLELCIEVVDNVQDAIDHIHKYGSSHTDIV TEDENTAEFFLQHVDVSACVFWNASTRFSDGYRFGGLGAEVGISTSRIHARGPVG LEGLLTTKWLLRGKDHVVSDFSEHGSLKYLHENLPIPQRNTN </p>
Research Area	Metabolism
Source	E.coli
Target Names	ALDH18A1
Expression Region	1-795aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

**Tag Info**


Tag-Free

Mol. Weight

87.3 kDa

Protein Length

Full Length

Image

(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.