





CD70 Recombinant Monoclonal Antibody

Product Code	CSB-RA004954MA1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P32970
Immunogen	Recombinant Human CD70 protein
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	Affinity-chromatography
Isotype	hIgG1
Clonality	Monoclonal
Product Type	Recombinant Antibody
Immunogen Species	Homo sapiens (Human)
Purity	Greater than 95% as determined by SEC-HPLC.
Target Names	CD70
Clone No.	5A4

Image



The Binding Activity of Human CD70 with Anti-CD70 recombinant antibody
Activity: Measured by its binding ability in a functional ELISA. Immobilized Human CD70 (CSB-MP004954HU1) at 2 µg/mL can bind Anti-CD70 recombinant antibody, the EC₅₀ is 1.642-2.005 ng/mL.



Greater than 95% as determined by SEC-HPLC.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.