





MUC17 Recombinant Monoclonal Antibody

Product Code	CSB-RA727848MA1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q685J3
Immunogen	Recombinant Human MUC17 protein
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	Affinity-chromatography
Isotype	Human IgG
Clonality	Monoclonal
Product Type	Recombinant Antibody
Immunogen Species	Homo sapiens (Human)
Purity	≥ 95% as determined by SEC-HPLC.
Research Area	Immunology
Target Names	MUC17
Clone No.	8D12

Image



The Binding Activity of Human MUC17 with Anti-MUC17 Recombinant Antibody
Activity: Measured by its binding ability in a functional ELISA. Immobilized Human MUC17 (CSB-MP727848HU) at 2 µg/mL can bind Anti-MUC17 recombinant antibody. The EC₅₀ is 0.9057-1.259 ng/mL.



The purity of MUC17 was greater than 95% as determined by SEC-HPLC

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.