





# MUC17 Recombinant Monoclonal Antibody

<b>Product Code</b>	CSB-RA727848MA1HU
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	Q685J3
<b>Immunogen</b>	Recombinant Human MUC17 protein
<b>Species Reactivity</b>	Human
<b>Tested Applications</b>	ELISA
<b>Form</b>	Liquid
<b>Conjugate</b>	Non-conjugated
<b>Storage Buffer</b>	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
<b>Purification Method</b>	Affinity-chromatography
<b>Isotype</b>	Human IgG
<b>Clonality</b>	Monoclonal
<b>Product Type</b>	Recombinant Antibody
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	Greater than 95% as determined by SEC-HPLC.
<b>Research Area</b>	Immunology
<b>Gene Names</b>	MUC17
<b>Clone No.</b>	8D12

## Image



The Binding Activity of Human MUC17 with Anti-MUC17 Recombinant Antibody  
Activity: Measured by its binding ability in a functional ELISA. Immobilized Human MUC17 (CSB-MP727848HU) at 2 µg/mL can bind Anti-MUC17 recombinant antibody. The EC<sub>50</sub> is 0.9057-1.259 ng/mL.



The purity of MUC17 was greater than 95% as determined by SEC-HPLC