




CD81 Recombinant Monoclonal Antibody

Product Code	CSB-RA004960MA1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P60033
Immunogen	Recombinant Human CD81 protein
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	Affinity-chromatography
Isotype	Human IgG1
Clonality	Monoclonal
Product Type	Recombinant Antibody
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology;Stem cells
Target Names	CD81
Clone No.	27H1

Image



The Binding Activity of Human CD81 with Anti-CD81 recombinant antibody
 Activity: Measured by its binding ability in a functional ELISA. Immobilized Human CD81(CSB-MP004960HUd7) at 2 µg/mL can bind Anti-CD81 recombinant antibody. the EC₅₀ is 4.166-5.578 ng/mL.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.