




ACVRL1 Recombinant Monoclonal Antibody

Product Code	CSB-RA001262MA1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P37023
Immunogen	Recombinant Human ACVRL1 protein
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	Affinity-chromatography
Isotype	Human IgG2
Clonality	Monoclonal
Product Type	Recombinant Antibody
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology;Stem cells
Gene Names	ACVRL1
Clone No.	6G7

Image



The Binding Activity of Human ACVRL1 with Anti-ACVRL1 recombinant antibody
Activity: Measured by its binding ability in a functional ELISA. Immobilized Human ACVRL1(CSB-BP001262HU1) at 2 µg/mL can bind Anti-ACVRL1 recombinant antibody, the EC₅₀ is 2.417-2.971ng/mL.