



CCR6 Recombinant Monoclonal Antibody

Product Code	CSB-RA004845MA2HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P51684
Immunogen	Recombinant Human CCR6 protein
Species Reactivity	Human
Tested Applications	Non-conjugated
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	Affinity-chromatography
Isotype	Human IgG1
Clonality	Monoclonal
Product Type	Recombinant Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	CCR6
Clone No.	7H5
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.