




Recombinant Mouse Carcinoembryonic antigen-related cell adhesion molecule 5 (Ceacam5)

Product Code	CSB-MP005165MO(A4)
Abbreviation	Recombinant Mouse Ceacam5 protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q3UKK2
Form	Liquid or Lyophilized powder
Storage Buffer	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose.
Product Type	Recombinant Mouse Carcinoembryonic antigen-related cell adhesion molecule 5 (Ceacam5)
Immunogen Species	Mus musculus (Mouse)
Sensitivity	Not Test
Purity	≥ 85% as determined by SDS-PAGE.
Sequence	QITIELEPPQVIEGENVLIRVNNLTENLITLAWFRGMRIKSPQIGQYTPATKVTVL GPGHSGRETLYSNGSLQIYNVTQEDIGFYSLRIINKHAEIVSITSIYLNVSSSLWT CEHPSPHAKLTIESVPPGISEGGSVLLL VKNLPQNLLSLFWYKGVIKVKFEIAR HIKATNSSVPGPAHTGRET VFSNGSLLLQEV MQSDTG FYTLRTMSTD LKDEVA HVQLYMDTYLLTCYHPLQVKIESLPQNVAVGKT VLLL VHNLPEDFQAFFWYKS AYRRD TYKIAEYKRAMDATILGSAYSSREFIYNNGSMLIIDVTEDDAGYF LLEILR EDLKIEKAYIQLHVNSFVSNSKDSASTARLSIESVPPSIVEGGSVLLL VHNLPKN LQSLFWYKGMIAEKKSELIQHIIATSSSLPGPAHSGRET VYNNNGSLLLQRVMQN DTQFYTLQTM DTLKYEVAHVQLQLDTS LSTWYHPLQVKIESLPRNVA VGKSV LFLVHNFPEVFRAFSWYKPAYKSHTSKIVEYHRFTDSATVGAAYRGIEVIFTNG SMVMIDVTEDDAGFYMLEILREDFKVEKAYVQLHVNSFVPNSKVS VSTARLSIE SVPPSIVGGESVLLL VHNLPKNLQSLFWYKGVIAEEKSELIQHIIATSSSLPGPA HSGRET VYSNGSLLLQRVMQN DTG FYTLT MSTD LKDEIAHVQLQLDTS TCCS LLTSDQLIIVPVPRNIAVGKSVLLLVCNVPKDVQTIFWYKSVYRTDIFKIAEYSRS MESTIWGLAHSGKEMVYTNGSLLIQNVTEHDAGLYMLEILHKDYKLERAHVQV HVNNPVSWPFVRVTDTTVRVQSSVFTCF SADPGVSIRWLFNKQSLQLTERM TLSPSKCQLSIDPVWREDAGKYQCEVSNPVSSKSSLPVRLAVIEE
Research Area	Cancer
Source	Mammalian cell
Target Names	Ceacam5
Expression Region	35-947aa



Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	C-terminal 10xHis-tagged
Mol. Weight	104 kDa
Protein Length	Full Length of Mature Protein

Image


(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
-----------------------	---

Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
-------------------	--