




KAAG1 Recombinant Monoclonal Antibody

Product Code	CSB-RA871385MA1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UBP8
Immunogen	Recombinant Human KAAG1 protein
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	Affinity-chromatography
Isotype	hIgG1
Clonality	Monoclonal
Product Type	Recombinant Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	KAAG1
Clone No.	2H3

Image



The Binding Activity of Human KAAG1 with Anti-KAAG1 recombinant antibody
Activity: Measured by its binding ability in a functional ELISA. Immobilized Human KAAG1(CSB-EP871385HU) at 2 µg/mL can bind Anti-KAAG1 recombinant antibody. The EC₅₀ is 2.040-2.284 ng/mL.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.