




CEACAM8 Recombinant Monoclonal Antibody

Product Code	CSB-RA005168MA1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P31997
Immunogen	Recombinant Human CEACAM8 protein
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	Affinity-chromatography
Isotype	hIgG1
Clonality	Monoclonal
Product Type	Recombinant Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	CEACAM8
Clone No.	22H9

Image



The Binding Activity of Human CEACAM8 with Anti-CEACAM8 Recombinant Antibody Activity: Measured by its binding ability in a functional ELISA. Immobilized Human CEACAM8 (CSB-MP005168HU) at 2 µg/mL can bind Anti-CEACAM8 recombinant antibody, the EC₅₀ is 19.38-21.68 ng/mL.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.