




IL6 Recombinant Monoclonal Antibody

Product Code	CSB-RA011664MA1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P05231
Immunogen	Recombinant Human IL6 protein
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	Affinity-chromatography
Isotype	hIgG1
Clonality	Monoclonal
Product Type	Recombinant Antibody
Immunogen Species	Homo sapiens (Human)
Gene Names	IL6
Clone No.	3C21

Image



The Binding Activity of Human IL6 with Anti-IL6 Recombinant Antibody
Activity: Measured by its binding ability in a functional ELISA. Immobilized Human IL6 (CSB-YP011664HU) at 2 µg/mL can bind Anti-IL6 recombinant antibody. The EC₅₀ is 35.80-41.82 ng/mL.