




DKK1 Recombinant Monoclonal Antibody

Product Code	CSB-RA006920MA1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O94907
Immunogen	Recombinant Human DKK1 protein
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	Affinity-chromatography
Isotype	hIgG1
Clonality	Monoclonal
Product Type	Recombinant Antibody
Immunogen Species	Homo sapiens (Human)
Gene Names	DKK1
Clone No.	27H11

Image



The Binding Activity of Human DKK1 with Anti-DKK1 Recombinant Antibody

Activity: Measured by its binding ability in a functional ELISA. Immobilized Human DKK1 (CSB-MP006920HU(A4)) at 2 µg/mL can bind Anti-DKK1 recombinant antibody, the EC₅₀ is 1.283-2.544 ng/mL.