




CEACAM5 Recombinant Monoclonal Antibody

Product Code	CSB-RA005165MA1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P06731
Immunogen	Recombinant Human CEACAM5 protein
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	Affinity-chromatography
Isotype	hIgG1
Clonality	Monoclonal
Product Type	Recombinant Antibody
Immunogen Species	Homo sapiens (Human)
Gene Names	CEACAM5
Clone No.	6G8

Image



The Binding Activity of Human CEACAM5 with Anti-CEACAM5 Recombinant Antibody Activity: Measured by its binding ability in a functional ELISA. Immobilized Human CEACAM5 (CSB-MP005165HU) at 2μg/mL can bind Anti-CEACAM5 recombinant antibody, the EC₅₀ is 0.8955-1.719 ng/mL.