




# flaB Antibody

|                            |  |
|----------------------------|--|
| <b>Product Code</b>        | CSB-PA527481XA01LOS  |
| <b>Storage</b>             | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.  |
| <b>Uniprot No.</b>         | O51941   |
| <b>Immunogen</b>           | Recombinant <i>Leptospira interrogans</i> serogroup Icterohaemorrhagiae serovar Lai flaB protein (1-283aa) |
| <b>Raised In</b>           | Rabbit   |
| <b>Species Reactivity</b>  | <i>Leptospira interrogans</i> serogroup Icterohaemorrhagiae serovar Lai                                    |
| <b>Tested Applications</b> | ELISA, WB; Recommended dilution: WB:1:500-1:5000   |
| <b>Form</b>                | Liquid   |
| <b>Conjugate</b>           | Non-conjugated   |
| <b>Storage Buffer</b>      | Preservative: 0.03% Proclin 300<br>Constituents: 50% Glycerol, 0.01M PBS, pH 7.4                           |
| <b>Purification Method</b> | >95%, Protein G purified   |
| <b>Isotype</b>             | IgG  |
| <b>Clonality</b>           | Polyclonal   |
| <b>Product Type</b>        | Polyclonal Antibody  |
| <b>Immunogen Species</b>   | <i>Leptospira interrogans</i> serogroup Icterohaemorrhagiae serovar Lai (strain 56601)                     |
| <b>Target Names</b>        | flaB   |

## Image



### Western Blot

Positive WB detected in: recombinant protein

All lanes: flaB Antibody at 1:1000

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 38.8 kDa

Observed band size: 38.8 kDa

## Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.