



Recombinant Human CCNA2&CDK2 Heterodimer Protein


Product Code	CSB-MP5050HU
Abbreviation	Recombinant Human CCNA2&CDK2 Heterodimer protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P20248/P24941
Form	Liquid or Lyophilized powder
Storage Buffer	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol.If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose.
Product Type	Recombinant Human CCNA2&CDK2 Heterodimer Protein
Immunogen Species	Homo sapiens (Human)
Sensitivity	Not Test
Purity	≥ 85% as determined by SDS-PAGE.
Sequence	MLGNSAPGPATREAGSALLALQQTALQEDQENINPEKAAPVQQPRTRAALAVL KSGNPRGLAQQQRPKTRRVAPLKDLVPNDEHVTVPWKANSKQPAFTIHVDE AEKEAQKKPAESQKIEREDALAFNSAISLPGPRKPLVPLDYPMDGSFESPHTM DMSIILEDEKPVSVNEVPDYHEDIHTYLREMEVKCKPKVGYMKKQPDITNSMR AILVDWLVEVGEEYKLQNETLHLAVNYIDRFLSSMSVLRGKLQLVGTAAMLLAS KFEEIYPPEVAEFVYITDDTYTKKQVLRMEHLVLKVLTFDLAAPTVMNQFLTQYFL HQQPANCKVESLAMFLGELSLIDADPYLKYLPSVIAGAAFHLALYTVTGSWPE SLIRKTGYTLESCLKPCLMDLHQTYLKAPQHAQQSIREKYKNSKYHGVSLNPPPE TLNL&MENFQKVEKIGEGTYGVVYKARNKLTGEVVALKKIRLDTETEGVPSTAI REISLLKELNHPNIVKLLDVIHTENKLYLVFEFLHQDLKKFMDASALTGIPLPLIKS YLFQLLQGLAFCHSHRVLHRDLKPQNLLINTEGAIKLADFGFLARAFGVPVRTYT HEVVTWYRAPEILLGCKYYSTAVDIWSLGCIFAEMVTRRALFPGDSEIDQLFRI FRTLGTDPDEVVWPGVTSMPPDYKPSFPKWARQDFSKVVPPLDEDGRSLLSQM LHYDPNKRISAKAALAHPPFFQDVTKPVPHLRL
Research Area	Cancer
Source	Mammalian cell
Target Names	CCNA2&CDK2
Expression Region	1-432aa&1-298aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	C-terminal mFc-tagged



Mol. Weight 77.9 kDa&35.3 kDa

Protein Length Heterodimer

Image



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.