




DXD Recombinant Monoclonal Antibody

Product Code	CSB-RA99697711m
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Immunogen	DXD-BSA conjugate
Species Reactivity	Chemical
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	Affinity-chromatography
Isotype	mIgG2a
Clonality	Monoclonal
Product Type	Recombinant Antibody
Research Area	Others
Target Names	DXD
Clone No.	8D12C11

Image



The Binding Activity of T-DXd(DS-8201) with Anti-DXD antibody
Activity: Measured by its binding ability in a functional ELISA. Immobilized T-DXd(DS-8201) at 2 µg/mL can bind Anti-DXD antibody, the EC₅₀ is 0.9828 to 1.446

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.