




DXD Recombinant Monoclonal Antibody

Product Code	CSB-RA99697711h
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Immunogen	DXD-BSA conjugate
Species Reactivity	Chemical
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	Affinity-chromatography
Isotype	Human IgG1
Clonality	Monoclonal
Product Type	Recombinant Antibody
Research Area	Others
Target Names	DXD
Clone No.	8D12C11

Image



The Binding Activity of DXD-BSA with Anti-DXD antibody

Activity: Measured by its binding ability in a functional ELISA. Immobilized DXD-BSA at 2 µg/mL can bind Anti-DXD antibody, the EC₅₀ is 22.68 to 27.25 ng/mL.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.