




SMO Antibody

Product Code	CSB-PA288930GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q99835
Immunogen	Recombinant Human Smoothened homolog protein (308-605AA)
Raised In	Rabbit
Species Reactivity	Human, Mouse, Rat
Tested Applications	ELISA, WB; Recommended dilution: WB:1:1000-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide and 50% glycerol pH 7.3
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	SMO

Image



Western Blot

Positive WB detected in: A549 whole cell lysate, MCF7 whole cell lysate, K562 whole cell lysate, 293T whole cell lysate, U251 whole cell lysate, HeLa whole cell lysate, Mouse Brain tissue lysate, Rat Brain tissue lysate

All lanes: BTRC antibody at 1:1000

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 69 kDa

Observed band size: 69 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.