




Recombinant Human Monocarboxylate transporter 6 (SLC16A5)-VLPs


Product Code	CSB-MP021412HU(A4)
Abbreviation	Recombinant Human SLC16A5 protein-VLPs
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	O15375
Form	Liquid or Lyophilized powder
Storage Buffer	Lyophilized from a 0.2 µm filtered PBS, 6% Trehalose, pH 7.4.
Product Type	Recombinant Human Monocarboxylate transporter 6 (SLC16A5)-VLPs
Immunogen Species	Homo sapiens (Human)
Sensitivity	Not Test
Purity	≥ 95% as determined by SEC-HPLC.
Sequence	MPQALERADGSWAWVLLATMVTQGLTLGFPTCIGIFFTELQWEFQASNSET SWFPSILTAVLHMAGPLCSILVGRFGCRVTVMLGGVLASLGMVASSFSHNLSQ LYFTAGFITGLGMCFSFQSSITVLGFYFVRRRVLANALASMGVSLGITLWPLLS RYLLENLGWRGTFLVFGGIFLHCCICGAIIRPVATSVAPETKECPPPPETPALG CLAACGRTIQRHLAFDILRHNTGYCVYILGVMWSVLGFPLPQVFLVPYAMWHS VDEQQAALLISIIGFSNIFLRPLAGLMAGRPAFASHRKYLFSLALLNGLTNLVCA ASGDFWVLVGYCLAYSVSMSGIGALIFQVLMDIVPMDQFPRALGLFTVLDGLA FLISPPLAGLLLDATNNFSYVFYMSSFFLISAALFMGGSFYALQKKEQGKQAVA ADALERDLFLEAKDGPQKQRSPEIMCQSSRQPRPAGVKNHLWGCPASSRTS HEWLLWPKAVLQAKQTALGWSNPT
Research Area	Cancer
Source	Mammalian cell
Target Names	SLC16A5
Expression Region	1-505aa
Notes	The VLPs are expressed from human 293 cells (HEK293). Mix the sample gently by repeatedly pipetting it up and down. Do not vortex. Repeated freezing and thawing is not recommended. Store the protein at -20°C/-80°C upon receiving it, and ensure to avoid repeated freezing and thawing, otherwise, it will affect the protein activity. The immunization strategy should be optimized (antigen dose, regimen and adjuvant).
Tag Info	C-terminal 10xHis-tagged(This tag can be tested only under denaturing conditions)
Mol. Weight	56.4 kDa


Protein Length

Full Length

Image


CSB-MP021412HU(A4) is detected by Mouse anti-6*His monoclonal antibody.



The purity of VLPs was greater than 95% as determined by SEC-HPLC

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. Aliquot for long-term storage at -80°C. Solubilize for 60 minutes at room temperature with occasional gentle mixing. Avoid vigorous shaking or vortexing.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.